



CHANDLER ENGINEERING

Model 4265HT

HIGH TEMPERATURE UCA

A Critical Tool for Oil Well Cement Testing

Testing the strength of standard cement cubes only provides information about cured strength. To accurately and effectively schedule a cementing job with a minimum WOC time, you need to determine how a slurry develops strength over time during curing. The Model 4265HT Ultrasonic Cement Analyzer (UCA) provides a determination of the strength development of a cement sample while it is being cured under very high downhole temperature and pressure conditions.

A Proven, Reliable Technique

Cement strength is determined by measuring the change in velocity of an ultrasonic signal transmitted through the cement specimen as it hardens. As the strength of the cement specimen increases, the ultrasonic signal's transit time through the sample decreases. Using proprietary empirical algorithms that have been proven throughout the Oil and Gas Industry, the relative strength is calculated.

Operational Simplicity

The cement slurry to be tested is prepared in accordance with API recommendations, then placed in the unit's temperature and pressure-controlled cell which simulates the curing conditions that are expected downhole. During testing, temperature is automatically controlled while pressure is manually set. For dynamic, precise pressure control, an automatic pressure system is available. All test data including the compressive strength results are then transferred to a computer, running Chandler Engineering Model 5270 Data Acquisition software.



FEATURES

- ✓ Real-Time Observation of Strength Development
- ✓ Used to predict WOC time
- ✓ Non-Destructive Method
- ✓ Uses Proven Industry-Standard Algorithms
- ✓ Chiller option for simulating low temperature cementing
- ✓ Chandler 5720 Data Acquisition Software
- ✓ Optional Automatic Pressure Controller
- ✓ Optional Model 4268ES Cement Expansion / Shrinkage cell

Model 4265HT

The software produces real-time graphs of the calculated strength, measured temperature and transit time. These graphs can be printed at any point during testing. All test data is recorded on the computer as the test progresses so that data is not lost in the event of a power failure. The Model 5270 software is capable of monitoring multiple instruments so only one computer is needed in the lab thereby saving valuable space.

Specifications

Maximum Operating Temperature
600°F (316°C)

Maximum Pressure
20,000 psi (138 MPa)

Utilities

Mains

208 – 240 VAC, 10A, 50/60 Hz, 1 Phase

Instrument Power

90 - 240 VAC, 1A, 50/60 Hz, 1 Phase

Water

20-80 psi / 140-550 kPa

Coolant

Water or Ethylene glycol solution

Air

Clean dry compressed air; 50 - 100 psi / 340 - 690 kPa

Drain

Suitable for hot water

Environmental

Indoor use, altitude up to 6562 ft (2000 m)

Ambient Temperature

60 - 122 °F (16 - 50 °C)

Storage Temperature

40 - 122 °F (5 - 50 °C)

Environment

Maximum Relatively Humidity

80% RH for temperatures up to 88 °F (31 °C)
decreasing linearly to 50% RH at 104 °F (40 °C)

Physical Dimensions (h x w x d)

21 in x 18 in x 19 in (53 cm x 46 cm x 48 cm)

Weight

130 lb (59 kg)

Shipping Information (h x w x d)

32 in x 26 in x 35 in (81 cm x 66 cm x 89 cm)

Weight

210 lb (95.5 kg)

Compliance

CE

Scan the below QR Code with your
phone to view product information
on our Website.



FOLLOW US ON:



AMETEK[®]**CHANDLER ENGINEERING**

2001 North Indianwood Avenue, Broken Arrow, OK 74012

Tel: +1 918-250-7200 • Fax: +1 918-459-0165

e-mail: chandler.sales@ametek.com • www.chandlereng.com

Houston Sales and Services

4903 W. Sam Houston Parkway, N., Suite A-400, Houston, TX 77041

Tel: +1 713-466-4900 • Fax: +1 713-849-1924